

Product specification

Blood line sealer MS 500M



Product information

Order number	87.7708.850
--------------	-------------

Product descriptions

Name	Blood line sealer MS 500M
Product number	87.7708.850
Product description	Mains operated tube sealer with hand-held sealer, RF cable (length: 1.8 m) and case

Product characteristics

Power consumption	200 W
Supply voltage	230 V

Size

Width of product	17 cm
Length of product	23.5 cm
Height of product	5.5 cm

Purity & certification

Accordance with MDD or IVDD	CE marked device
CE certified	CE
Plasticiser	yes
Lot or serial number	no

Packaging

Piece(s) / pallet	84
Height of case	125 mm
Width of case	358 mm
Depth of case	388 mm
Volume of case	0.0174 cbm

Weight of case

3.8 kg

EAN inner box

4038917361637

EAN case

4038917361620

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

