



Product specification

Micro sample tube Fluoride/heparin FH, 1.3 ml, screw cap, ISO



Product information

Product number	41.1394.105
----------------	-------------

Product descriptions

Name	Micro sample tube Fluoride/heparin FH, 1.3 ml, screw cap, ISO
Product number	41.1394.105
Product description	Micro sample tube Fluoride/heparin FH, venous blood collection, preparation: Fluoride + heparin, nominal volume: 1.3 ml, (LxØ) with cap: 47 x 10.8 mm, screw cap, cap: grey, colour code ISO, with paper label, label/print: grey, skirted conical base, 100 piece(s)/bag

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	Polypropylene (PP)
Color of cap	grey

Size

Diameter	10.8 mm
Length excluding cap	45 mm

Length including cap	47 mm
Length of product	47 mm
Sample volume	1.3 ml

Product characteristics

Base shape	skirted conical base
Label/Print	With paper label
Color of print/label	grey
Graduation	yes
Preparation concentration	1.0 mg/ml blood + 16 IU/ml blood
Closure type	Screw cap
Type of collection	venous
Centrifugation recommended	2000 x g
Centrifugation time	5 Min

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	96000
Height of case	140 mm
Width of case	278 mm
Depth of case	348 mm
Volume of case	0.0135 cbm
Weight of case	1.91 kg
EAN inner box	4038917128353
EAN case	4038917021678

This document was prepared by EDP support
and is valid without signature.



SARSTEDT