



# Product specification

## Rack, styrofoam, format: 5 x 4



### Product information

Product number	95.64.259
----------------	-----------

### Product descriptions

Name	Rack, styrofoam, format: 5 x 4
Product number	95.64.259
Product description	Rack, material: styrofoam, white, format: 5 x 4, for 20 collection tubes, 8 piece(s)/case

### Material & colours

Color of product	white
------------------	-------

### Product characteristics

Storage capacity	20
------------------	----

### Purity & certification

Product category	no medical device   IVD
Batched	yes

### Packaging

Piece(s) / case	56
Piece(s) / pallet	1792
Volume of case	0.0559 cbm
Weight of case	1.36 kg

EAN inner box

4038917386883

EAN case

4038917386876

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT