



Product specification

Silicone strip, length: 240 mm, translucent



Product information

Product number	93.1647.050
----------------	-------------

Product descriptions

Name	Silicone strip, length: 240 mm, translucent
Product number	93.1647.050
Product description	Silicone strip, for DishRack 50, length: 240 mm, material: silicone, translucent, 1 piece(s)/bag

Material & colours

Color of product	translucent
------------------	-------------

Size

Width of product	65 mm
Length of product	240 mm

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Piece(s) / case	30
-----------------	----

Piece(s) / pallet	4860
Height of inner box	27 mm
Width of inner box	65 mm
Depth of inner box	240 mm
Height of case	190 mm
Width of case	193 mm
Depth of case	253 mm
Volume of case	0.0093 cbm
Weight of case	2 kg
EAN inner box	4038917132640
EAN case	4038917040648

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT