



## Product specification

**DishRack, height: 240 mm, green, for 52 petri dishes with 92 mm Ø**



## Product information

Product number	93.1647.004
----------------	-------------

## Product descriptions

Name	DishRack, height: 240 mm, green, for 52 petri dishes with 92 mm Ø
Product number	93.1647.004
Product description	DishRack, height: 240 mm, material: PC, green, for 52 petri dishes with 92 mm Ø, 1 piece(s)/case

## Material & colours

Color of product	green
------------------	-------

## Size

Height of product	240 mm
-------------------	--------

## Product characteristics

Storage capacity	52
------------------	----

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / pallet	162
Height of inner box	190 mm
Width of inner box	193 mm
Depth of inner box	253 mm
Height of case	190 mm
Width of case	193 mm
Depth of case	253 mm
Volume of case	0.0093 cbm
Weight of case	0.8 kg
EAN inner box	4038917132633
EAN case	4038917030540

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT