



## Product specification

**DishRack, height: 370 mm, green, for 88 petri dishes with 92 mm Ø**



## Product information

Product number	93.1646.004
----------------	-------------

## Product descriptions

Name	DishRack, height: 370 mm, green, for 88 petri dishes with 92 mm Ø
Product number	93.1646.004
Product description	DishRack, height: 370 mm, material: PC, green, for 88 petri dishes with 92 mm Ø, 1 piece(s)/case

## Material & colours

Color of product	green
------------------	-------

## Size

Height of product	370 mm
-------------------	--------

## Product characteristics

Storage capacity	88
------------------	----

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / pallet	108
Height of inner box	190 mm
Width of inner box	193 mm
Depth of inner box	383 mm
Height of case	190 mm
Width of case	193 mm
Depth of case	383 mm
Volume of case	0.014 cbm
Weight of case	1.35 kg
EAN inner box	4038917132589
EAN case	4038917030496

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT