




Product specification

Inoculation spreader, PS, blue, sterile



Product information

Product number	86.1569.005
----------------	-------------

Product descriptions

Name	Inoculation spreader, PS, blue, sterile
Product number	86.1569.005
Product description	Inoculation spreader, material: PS, blue, for spreading a bacterial suspension, sterile, 4 piece(s)/blister

Material & colours

Product material	Polystyrene (PS)
Color of product	blue

Size

Width of product	40 mm
Length of product	185 mm

Product characteristics

Example of use	for spreading a bacterial suspension
----------------	--------------------------------------

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation

Batched

yes

Packaging

Minimum order quantity	500
Type of inner packaging	blister
Piece(s) / case	500
Piece(s) / pallet	33000
Height of case	155 mm
Width of case	268 mm
Depth of case	438 mm
Volume of case	0.0182 cbm
Weight of case	1.7 kg
EAN inner box	4038917130998
EAN case	4038917023566

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT