



Product specification

Storage box, slip-on lid, cardboard, format: 8 x 8, for 64 collection tubes



Product information

Product number	95.64.964
----------------	-----------

Product descriptions

Name	Storage box, Slip-on lid, Cardboard, Format: 8 x 8, for 64 Collection tubes
Product number	95.64.964
Product description	Storage box, slip-on lid, material: cardboard, white, format: 8 x 8, for 64 collection tubes, suitable for tubes Ø 15 mm, 1 piece(s)/case

Material & colours

Color of product	white
Closure material	Cardboard
Color of cap	white

Size

Width of product	135 mm
Length of product	135 mm
Height of product	45 mm

Product characteristics

Label/Print	With print
-------------	------------

Storage capacity	64
------------------	----

Closure type	Slip-on lid
--------------	-------------

Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
--------------	------------

Batched	yes
---------	-----

Packaging

Piece(s) / case	120
-----------------	-----

Piece(s) / pallet	1920
-------------------	------

Volume of case	0.1134 cbm
----------------	------------

Weight of case	10.95 kg
----------------	----------

EAN inner box	4038917247412
---------------	---------------

EAN case	4038917247405
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

