




## Product specification

### Inoculation needle, PS, orange, sterile



## Product information

Product number	86.1568.050
----------------	-------------

## Product descriptions

Name	Inoculation needle, PS, orange, sterile
Product number	86.1568.050
Product description	Inoculation needle, material: PS, orange, sterile, 50 piece(s)/blister

## Material & colours

Color of product	orange
------------------	--------

## Size

Width of product	4 mm
Length of product	197 mm

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

## Packaging

Type of inner packaging	blister
-------------------------	---------

Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	264000
Height of inner box	101 mm
Width of inner box	161 mm
Depth of inner box	216 mm
Height of case	120 mm
Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0092 cbm
Weight of case	2.79 kg
EAN inner box	4038917130974
EAN case	4038917023542

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT