




Product specification

**PCR plate without skirt, 96 well, white, High Profile, 200 µl,
PCR Performance Tested, PP**



Product information

Product number	72.1978.010
----------------	-------------

Product descriptions

Name	PCR plate without skirt, 96 well, white, High Profile, 200 µl, PCR Performance Tested, PP
Product number	72.1978.010
Product description	PCR plate without skirt, 96 well, white, High Profile, 200 µl, PCR Performance Tested, material: PP, optimised for qPCR, 10 piece(s)/minigrip bag

Material & colours

Color of product	white
Profile	High Profile

Size

Diameter	5.51 mm
Width of product	76 mm
Length of product	119.9 mm
Height of product	21 mm
Max. volume	200 µl
Sample volume	0.2 ml

Product characteristics

Number of wells	96
Construction of the edges of the plate	without skirt
Base shape	conical
Closure type	Film
Autoclavable (max. °C)	121 °C
Suitable for qPCR	yes
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	100
Piece(s) / pallet	5600
Height of inner box	131 mm
Width of inner box	230 mm
Depth of inner box	270 mm
Height of case	247 mm
Width of case	280 mm
Depth of case	282 mm
Volume of case	0.0195 cbm
Weight of case	2.27 kg
EAN inner box	4038917557665
EAN case	4038917557658

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT