




Product specification

PCR plate half skirt, 96 well, transparent, High Profile, 200 µl, Low DNA-binding, PCR Performance T



Product information

Product number	72.1979.700
----------------	-------------

Product descriptions

Name	PCR plate half skirt, 96 well, transparent, High Profile, 200 µl, Low DNA-binding, PCR Performance Tested, PP
Product number	72.1979.700
Product description	PCR plate half skirt, 96 well, transparent, High Profile, 200 µl, Low DNA-binding, PCR Performance Tested, material: PP, 10 piece(s)/minigrip bag

Material & colours

Color of product	transparent
Profile	High Profile

Size

Diameter	4.95 mm
Width of product	86 mm
Length of product	126.1 mm
Height of product	19.33 mm
Max. volume	200 µl

Sample volume	0.2 ml
---------------	--------

Product characteristics

Number of wells	96
Construction of the edges of the plate	half skirt
Base shape	conical
Suitable for qPCR	yes
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	PCR Performance Tested
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	100
Piece(s) / pallet	5600
Height of inner box	131 mm
Width of inner box	230 mm
Depth of inner box	270 mm
Height of case	247 mm
Width of case	280 mm
Depth of case	282 mm
Volume of case	0.0195 cbm
Weight of case	2.29 kg
EAN inner box	4038917574631
EAN case	4038917574624

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT