


Product specification

PCR plate full skirt, 96 well, transparent, Low Profile, 100 μ l, Low DNA-binding, PCR Performance Te



Product information

Order number	72.1980.700
--------------	-------------

Product descriptions

Name	PCR plate full skirt, 96 well, transparent, Low Profile, 100 μ l, Low DNA-binding, PCR Performance Tested, PP
Product number	72.1980.700
Product description	PCR plate full skirt, 96 well, transparent, Low Profile, 100 μ l, Low DNA-binding, PCR Performance Tested, material: PP, 10 piece(s)/bag

Product characteristics

Number of wells	96
Construction of the edges of the plate	full skirt
Base shape	conical
Suitable for qPCR	yes
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22 $^{\circ}$ C
Centrifugation time	5 Min

Size

Diameter	5.5 mm
Width of product	85.5 mm
Length of product	127.8 mm
Height of product	15.5 mm
Max. volume	100 μ l
Sample volume	0.1 ml

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Profile	Low Profile

Purity & certification

Accordance with MDD or IVDD	no medical device IVD
Purity standard	PCR Performance Tested
Sterilisation	Electron irradiation
Lot or serial number	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	100
Piece(s) / pallet	5600
Height of case	247 mm
Width of case	280 mm
Depth of case	282 mm
Volume of case	0.0195 cbm
Weight of case	2.34 kg
EAN inner box	4038917574679
EAN case	4038917574662

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

