




Product specification

PCR plate full skirt, 96 well, transparent, Low Profile, 100 µl, Low DNA-binding, PCR Performance Te



Product information

Product number	72.1980.700
----------------	-------------

Product descriptions

Name	PCR plate full skirt, 96 well, transparent, Low Profile, 100 µl, Low DNA-binding, PCR Performance Tested, PP
Product number	72.1980.700
Product description	PCR plate full skirt, 96 well, transparent, Low Profile, 100 µl, Low DNA-binding, PCR Performance Tested, material: PP, 10 piece(s)/minigrip bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Profile	Low Profile

Size

Diameter	5.5 mm
Width of product	85.5 mm
Length of product	127.8 mm
Height of product	15.5 mm

Max. volume	100 µl
Sample volume	0.1 ml

Product characteristics

Number of wells	96
Construction of the edges of the plate	full skirt
Base shape	conical
Suitable for qPCR	yes
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	100
Piece(s) / inner box	50
Piece(s) / case	100
Piece(s) / pallet	5600
Height of inner box	131 mm
Width of inner box	230 mm
Depth of inner box	270 mm
Height of case	247 mm
Width of case	280 mm
Depth of case	282 mm
Volume of case	0.0195 cbm
Weight of case	2.34 kg
EAN inner box	4038917574679
EAN case	4038917574662

This document was prepared by EDP support
and is valid without signature.



SARSTEDT