

Product specification

Storage box, slip-on lid, cardboard, format: 7 x 7, for 49 collection tubes



Product information

Product number	95.64.949
----------------	-----------

Product descriptions

Name	Storage box, Slip-on lid, Cardboard, Format: 7 x 7, for 49 Collection tubes
------	---

Product number	95.64.949
----------------	-----------

Product description	Storage box, slip-on lid, material: cardboard, white, format: 7 x 7, for 49 collection tubes, suitable for tubes up to 17 mm Ø (15 ml tubes), 1 piece(s)/case
---------------------	---

Material & colours

Color of product	white
------------------	-------

Closure material	Cardboard
------------------	-----------

Color of cap	white
--------------	-------

Size

Width of product	135 mm
------------------	--------

Length of product	135 mm
-------------------	--------

Height of product	105 mm
-------------------	--------

Product characteristics

Storage capacity	49
Closure type	Slip-on lid

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / case	60
Piece(s) / pallet	720
Height of inner box	105 mm
Width of inner box	135 mm
Depth of inner box	135 mm
Height of case	560 mm
Width of case	440 mm
Depth of case	560 mm
Volume of case	0.1132 cbm
Weight of case	8.35 kg
EAN inner box	4038917095792
EAN case	4038917561686

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

