




Product specification

Inoculation loop, 10 µl, PS, blue, sterile



Product information

Product number	86.1562.050
----------------	-------------

Product descriptions

Name	Inoculation loop, 10 µl, PS, blue, sterile
Product number	86.1562.050
Product description	Inoculation loop, volume: 10 µl, material: PS, blue, sterile, 48 piece(s)/blister

Material & colours

Product material	Polystyrene (PS)
Color of product	blue

Size

Width of product	6 mm
Length of product	200 mm

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order quantity	1920
Type of inner packaging	blister
Piece(s) / inner box	960
Piece(s) / case	1920
Piece(s) / pallet	253440
Height of inner box	101 mm
Width of inner box	161 mm
Depth of inner box	216 mm
Height of case	120 mm
Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0092 cbm
Weight of case	2.5805 kg
EAN inner box	4038917130882
EAN case	4038917073127

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT