




Product specification

Serological pipette, without tip, plugged, 10 ml, sterile, 10 piece(s)/bag



Product information

Product number	86.1688.010
----------------	-------------

Product descriptions

Name	Serological pipette, without tip, plugged, 10 ml, sterile, 10 piece(s)/Bag
Product number	86.1688.010
Product description	Serological pipette, without tip, plugged, 10 ml, material: PS, sterile, 10 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Diameter	10 mm
Length of product	335 mm

Product characteristics

Design	without tip, plugged
Graduation	yes

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	10
Piece(s) / case	500
Piece(s) / pallet	27000
Height of case	270 mm
Width of case	250 mm
Depth of case	413 mm
Volume of case	0.0279 cbm
Weight of case	4.97 kg
EAN inner box	4038917133081
EAN case	4038917024693

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT