



Product specification

Microvette® APT 250 EDTA K2E, 250 µl, cap violet, cap, round base



Product information

Product number	20.1331.100
----------------	-------------

Product descriptions

Name	Microvette® APT 250 EDTA K2E, 250 µl, cap violet, cap, round base
Product number	20.1331.100
Product description	Microvette® APT 250 EDTA K2E, capillary blood collection, preparation: K2 EDTA, K2E, nominal volume: 250 µl, Ø: 13 mm, with plastic label, cap violet, cap, with pierceable membrane, round base, inner receptacle: cylindrical, incl. quick-release cap, non-sterile, 50 piece(s)/case

Material & colours

Color of product	transparent
Closure material	Polyethylene (PE)
Color of cap	violet

Size

Diameter	13 mm
Length excluding cap	75 mm
Sample volume	250 µl

Product characteristics

Base shape	Round base
Label/Print	With plastic label
Color of print/label	transparent/violet
Shape of inner tube	cylindrical
Cap	Membrane cap
Closure type	Cap
Type of collection	capillary
Pierceable	yes

Purity & certification

Satisfies the requirement	P 650, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	non-sterile
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	32000
Height of inner box	72 mm
Width of inner box	128 mm
Depth of inner box	201 mm
Height of case	221 mm
Width of case	268 mm
Depth of case	378 mm
Volume of case	0.0224 cbm
Weight of case	4.01 kg
EAN inner box	4038917543026
EAN case	4038917543019

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT