




## Product specification

**Pipette tip, 200 µl, transparent, Biosphere® plus, 96 piece(s)/SingleRefill**



## Product information

|                |             |
|----------------|-------------|
| Product number | 70.3030.305 |
|----------------|-------------|

## Product descriptions

|                     |  |
|---------------------|--|
| Name                | Pipette tip, 200 µl, transparent, Biosphere® plus, 96 piece(s)/SingleRefill  |
| Product number      | 70.3030.305  |
| Product description | Pipette tip, working volume: 200 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/SingleRefill |

## Material & colours

|                  |             |
|------------------|-------------|
| Color of product | transparent |
|------------------|-------------|

## Size

|                   |        |
|-------------------|--------|
| Length of product | 51 mm  |
| Max. volume       | 200 µl |

## Product characteristics

|        |    |
|--------|----|
| Filter | no |
|--------|----|

|                        |        |
|------------------------|--------|
| Graduation             | yes    |
| Autoclavable (max. °C) | 121 °C |

## Purity & certification

|                  |                                      |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD)            |
| CE certified     | CE - manufacturer's self-declaration |
| Purity standard  | Biosphere plus                       |
| Sterilisation    | Gas sterilisation                    |
| Batched          | yes                                  |

## Packaging

|                      |               |
|----------------------|---------------|
| Piece(s) / inner box | 480           |
| Piece(s) / case      | 1920          |
| Piece(s) / pallet    | 115200        |
| Height of inner box  | 148 mm        |
| Width of inner box   | 105 mm        |
| Depth of inner box   | 291 mm        |
| Height of case       | 170 mm        |
| Width of case        | 307 mm        |
| Depth of case        | 437 mm        |
| Volume of case       | 0.0228 cbm    |
| Weight of case       | 2.566 kg      |
| EAN inner box        | 4038917847032 |
| EAN case             | 4038917847049 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT