




## Product specification

**Serological pipette, without tip, plugged, 5 ml, sterile, 10 piece(s)/bag**



### Product information

Product number	86.1687.010
----------------	-------------

### Product descriptions

Name	Serological pipette, without tip, plugged, 5 ml, sterile, 10 piece(s)/Bag
Product number	86.1687.010
Product description	Serological pipette, without tip, plugged, 5 ml, material: PS, sterile, 10 piece(s)/bag

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Diameter	8 mm
Length of product	340 mm
Sample volume	5 ml

### Product characteristics

Design	without tip, plugged
Label/Print	with printed graduations
Color of print/label	black
Graduation	yes

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

## Packaging

Piece(s) / inner box	10
Piece(s) / case	500
Piece(s) / pallet	40500
Height of case	195 mm
Width of case	250 mm
Depth of case	413 mm
Volume of case	0.0201 cbm
Weight of case	3.88 kg
EAN inner box	4038917133074
EAN case	4038917024686

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT