




## Product specification

**Demeter pipette, without tip, plugged, 1.1 ml, sterile, 25 piece(s)/bag**



## Product information

Product number	86.1686.225
----------------	-------------

## Product descriptions

Name	Demeter pipette, without tip, plugged, 1.1 ml, sterile, 25 piece(s)/Bag
Product number	86.1686.225
Product description	Demeter pipette, without tip, plugged, 1.1 ml, material: PS, sterile, 25 piece(s)/bag

## Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

## Size

Diameter	4 mm
Length of product	270 mm

## Product characteristics

Design	without tip, plugged
Graduation	yes

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	25
Piece(s) / case	1000
Piece(s) / pallet	108000
Height of case	150 mm
Width of case	250 mm
Depth of case	413 mm
Volume of case	0.0155 cbm
Weight of case	4.04 kg
EAN inner box	4038917133050
EAN case	4038917024679

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT