




## Product specification

**Pipette tip, 250 µl, transparent, Biosphere® plus, 96 piece(s)/box**



## Product information

Product number	70.3031.205
----------------	-------------

## Product descriptions

Name	Pipette tip, 250 µl, transparent, Biosphere® plus, 96 piece(s)/Box
Product number	70.3031.205
Product description	Pipette tip, working volume: 250 µl, transparent, filling level rings, Biosphere® plus, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Length of product	51 mm
Max. volume	250 µl

## Product characteristics

Filter	no
Graduation	yes

Autoclavable (max. °C)	121 °C
------------------------	--------

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Batched	yes

## Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	96000
Height of inner box	166 mm
Width of inner box	129 mm
Depth of inner box	348 mm
Height of case	245 mm
Width of case	325 mm
Depth of case	348 mm
Volume of case	0.0277 cbm
Weight of case	3.04 kg
EAN inner box	4038917845793
EAN case	4038917845809

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT