




## Product specification

**Serological pipette, with tip, plugged, 25 ml, sterile, 20 piece(s)/bag**



### Product information

Product number	86.1685.020
----------------	-------------

### Product descriptions

Name	Serological pipette, with tip, plugged, 25 ml, sterile, 20 piece(s)/bag
Product number	86.1685.020
Product description	Serological pipette, with tip, plugged, 25 ml, coding colour: red, material: PS, sterile, 20 piece(s)/bag

### Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

### Size

Diameter	15 mm
Length of product	345 mm

### Product characteristics

Design	with tip, plugged
Color of print/label	red
Graduation	yes

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

## Packaging

Minimum order quantity	200
Piece(s) / inner box	20
Piece(s) / case	200
Piece(s) / pallet	10800
Height of case	270 mm
Width of case	250 mm
Depth of case	413 mm
Volume of case	0.0279 cbm
Weight of case	3.61 kg
EAN inner box	4038917133029
EAN case	4038917040785

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT