



## Product specification

### Screw cap micro tube, 2 ml, sterile



## Product information

Product number	72.694.325
----------------	------------

## Product descriptions

Name	Screw cap micro tube, 2 ml, sterile
Product number	72.694.325
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, cap: green, cap assembled, with printed writing space, with graduations, sterile, 100 piece(s)/bag

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
O-Ring material	Ethylene propylene diene rubber (EPDM)
Closure material	Polypropylene (PP)
Color of cap	green

## Size

Diameter	10.8 mm
Length of product	44 mm

## Product characteristics

Base shape	skirted conical base
------------	----------------------

Label/Print	With print
Color of print/label	white
Graduation	yes
Knurls	yes
Centrifugation max (RCF)	20000 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	108000
Height of case	138 mm
Width of case	280 mm
Depth of case	351 mm
Volume of case	0.0136 cbm
Weight of case	1.82 kg
EAN inner box	4038917926614
EAN case	4038917926621

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT