




## Product specification

**Serological pipette, with tip, plugged, 10 ml, sterile, 25 piece(s)/bag**



### Product information

Product number	86.1254.025
----------------	-------------

### Product descriptions

Name	Serological pipette, with tip, plugged, 10 ml, sterile, 25 piece(s)/bag
Product number	86.1254.025
Product description	Serological pipette, with tip, plugged, 10 ml, coding colour: orange, material: PS, sterile, 25 piece(s)/bag

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Diameter	10 mm
Length of product	335 mm
Sample volume	10 ml

### Product characteristics

Design	with tip, plugged
Label/Print	with printed graduations
Color of print/label	orange/black
Graduation	yes

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

## Packaging

Piece(s) / inner box	25
Piece(s) / case	500
Piece(s) / pallet	27000
Height of case	270 mm
Width of case	250 mm
Depth of case	413 mm
Volume of case	0.0279 cbm
Weight of case	5.26 kg
EAN inner box	4038917127318
EAN case	4038917021135

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT