



## Product specification

### Screw cap, natural, sterile, suitable for screw cap micro tubes



## Product information

Product number	65.716.320
----------------	------------

## Product descriptions

Name	Screw cap, natural, sterile, suitable for Screw cap micro tubes
Product number	65.716.320
Product description	Screw cap, natural, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag

## Material & colours

Color of product	natural
O-Ring material	Ethylene propylene diene rubber (EPDM)

## Product characteristics

Cap	flat cap
Closure type	Screw cap

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE
Purity standard	sterile

Sterilisation

Electron irradiation

Batched

yes

## Packaging

Piece(s) / inner box

500

Piece(s) / case

5000

Piece(s) / pallet

540000

Height of case

138 mm

Width of case

280 mm

Depth of case

351 mm

Volume of case

0.0136 cbm

Weight of case

2.05 kg

EAN inner box

4038917923323

EAN case

4038917923330

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT