



Product specification

Screw cap, red, sterile, suitable for screw cap micro tubes



Product information

Product number	65.716.323
----------------	------------

Product descriptions

Name	Screw cap, red, sterile, suitable for Screw cap micro tubes
Product number	65.716.323
Product description	Screw cap, red, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag

Material & colours

Color of product	red
O-Ring material	Ethylene propylene diene rubber (EPDM)

Product characteristics

Cap	flat cap
Closure type	Screw cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation

Batched	yes
---------	-----

Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	540000
Height of case	138 mm
Width of case	280 mm
Depth of case	351 mm
Volume of case	0.0136 cbm
Weight of case	2.05 kg
EAN inner box	4038917923415
EAN case	4038917923422

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT