



## Product specification

### Screw cap, orange, sterile, suitable for screw cap micro tubes



#### Product information

Product number	65.716.326
----------------	------------

#### Product descriptions

Name	Screw cap, orange, sterile, suitable for Screw cap micro tubes
Product number	65.716.326
Product description	Screw cap, orange, sterile, suitable for screw cap micro tubes, 500 piece(s)/double bag

#### Material & colours

Product material	Polypropylene (PP)
Color of product	orange
O-Ring material	Ethylene propylene diene rubber (EPDM)

#### Product characteristics

Cap	flat cap
Closure type	Screw cap

#### Purity & certification

Product category	diagnostic device, CE
CE certified	CE

Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	5000
Piece(s) / pallet	540000
Height of case	138 mm
Width of case	280 mm
Depth of case	351 mm
Volume of case	0.0136 cbm
Weight of case	2.05 kg
EAN inner box	4038917923507
EAN case	4038917923514

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT