



Product specification

Screw cap tube, 5 ml, round base, transparent



Product information

Product number	63.504.027
----------------	------------

Product descriptions

Name	Screw cap tube, 5 ml, round base, transparent
Product number	63.504.027
Product description	Screw cap tube, working volume: 5 ml, (LxØ): 75 x 13 mm, round base, PP, transparent, cap assembled, with plastic label, label/print: transparent/green, "CSF/Liquor", with graduations, cap assembled, natural, sterile, 1 piece(s)/blister

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

Size

Diameter	13 mm
Length excluding cap	75 mm
Length including cap	83 mm
Length of product	75 mm

Product characteristics

Base shape	Round base
Label/Print	With plastic label
Color of print/label	transparent/green
Graduation	yes
Closure type	Screw cap
Example of use	CSF, collection, transport and storage of CSF/liquor

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes
Safety data sheet	no

Packaging

Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	13500
Height of inner box	171 mm
Width of inner box	151 mm
Depth of inner box	351 mm
Height of case	195 mm
Width of case	368 mm
Depth of case	778 mm
Volume of case	0.0558 cbm
Weight of case	3.56 kg
EAN inner box	4038917408561
EAN case	4038917408554

This document was prepared by EDP support
and is valid without signature.



SARSTEDT