



Product specification

Cell culture plate, 96 well, surface: BIOFLOAT™, round base



Product information

Product number	83.3925.400
----------------	-------------

Product descriptions

Name	Cell culture plate, 96 well, surface: BIOFLOAT™, round base
Product number	83.3925.400
Product description	Cell culture plate, 96 well, material: PS, surface: BIOFLOAT™, for spheroid culture, round base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	Polystyrene (PS)

Size

Width of product	85.2 mm
Length of product	127.8 mm
Height of product	16.55 mm
Max. volume	0.31 ml

Number of wells	96
Base shape	Round base
Vessel type	Plate
Closure type	Slip-on lid
Cell type	spheroid culture

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / case	24
Piece(s) / pallet	1440
Height of inner box	80 mm
Width of inner box	140 mm
Depth of inner box	205 mm
Height of case	255 mm
Width of case	257 mm
Depth of case	322 mm
Volume of case	0.0211 cbm
Weight of case	2.357 kg
EAN inner box	4038917923965
EAN case	4038917923972

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT