




## Product specification

**Pipette tip, 10 µl, transparent, PCR Performance Tested,  
96 piece(s)/box**



## Product information

Product number	70.3010.200
----------------	-------------

## Product descriptions

Name	Pipette tip, 10 µl, transparent, PCR Performance Tested, 96 piece(s)/box
Product number	70.3010.200
Product description	Pipette tip, working volume: 10 µl, transparent, filling level rings, PCR Performance Tested, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/box

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Length of product	31 mm
Max. volume	10 µl

## Product characteristics

Filter	no
Graduation	yes

Autoclavable (max. °C)	121 °C
------------------------	--------

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	96000
Height of inner box	148 mm
Width of inner box	111 mm
Depth of inner box	330 mm
Height of case	245 mm
Width of case	325 mm
Depth of case	348 mm
Volume of case	0.0277 cbm
Weight of case	2.74 kg
EAN inner box	4038917894708
EAN case	4038917894715

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT