




Product specification

Aspiration pipette, with tip, unplugged, 2 ml, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic,



Product information

Product number	86.1252.011
----------------	-------------

Product descriptions

Name	Aspiration pipette, with tip, unplugged, 2 ml, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister
Product number	86.1252.011
Product description	Aspiration pipette, with tip, unplugged, 2 ml, material: PS, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Diameter	6 mm
Length of product	270 mm

Product characteristics

Design	with tip, unplugged
Color of print/label	green

Graduation	no
------------	----

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	54000
Height of case	290 mm
Width of case	250 mm
Depth of case	420 mm
Volume of case	0.0323 cbm
Weight of case	5.03 kg
EAN inner box	4038917127202
EAN case	4038917043311

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT