




Product specification

Filter tip, 2.5 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/box



Product information

Product number	70.3010.285
----------------	-------------

Product descriptions

Name	Filter tip, 2.5 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/Box
Product number	70.3010.285
Product description	Filter tip, working volume: 2.5 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/box

Material & colours

Color of product	transparent
Filter material	Polyethylene (PE)

Size

Length of product	31 mm
Max. volume	2.5 µl

Filter	yes
Graduation	yes

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	96000
Height of inner box	148 mm
Width of inner box	111 mm
Depth of inner box	330 mm
Height of case	245 mm
Width of case	325 mm
Depth of case	348 mm
Volume of case	0.0277 cbm
Weight of case	2.7405 kg
EAN inner box	4038917894661
EAN case	4038917894678

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT