




## Product specification

**Filter tip, 2.5 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/SingleRefill**



## Product information

Product number	70.3010.385
----------------	-------------

## Product descriptions

Name	Filter tip, 2.5 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/SingleRefill
Product number	70.3010.385
Product description	Filter tip, working volume: 2.5 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/SingleRefill

## Material & colours

Color of product	transparent
Filter material	Polyethylene (PE)

## Size

Length of product	31 mm
Max. volume	2.5 µl

Filter	yes
Graduation	yes

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Batched	yes

## Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	115200
Height of inner box	148 mm
Width of inner box	105 mm
Depth of inner box	291 mm
Height of case	170 mm
Width of case	307 mm
Depth of case	437 mm
Volume of case	0.0228 cbm
Weight of case	2.16 kg
EAN inner box	4038917894913
EAN case	4038917894920

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT