


Product specification

Serological pipette, plugged, 2 ml, sterile, 25 piece(s)/bag



Product information

Order number	86.1252.025
--------------	-------------

Product descriptions

Name	Serological pipette, plugged, 2 ml, sterile, 25 piece(s)/Bag
Product number	86.1252.025
Product description	Serological pipette, plugged, 2 ml, coding colour: green, material: PS, sterile, 25 piece(s)/bag

Product characteristics

Design	plugged
Color of print/label	green
Graduation	yes

Size

Diameter	6 mm
Length of product	270 mm

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Purity & certification

Accordance with MDD or IVDD	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Lot or serial number	yes

Piece(s) / inner box	25
Piece(s) / case	1000
Piece(s) / pallet	81000
Volume of case	0.0194 cbm
Weight of case	4.21 kg
EAN inner box	4038917127219
EAN case	4038917021098

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

