




Product specification

Serological pipette, with tip, plugged, 2 ml, sterile, 25 piece(s)/bag



Product information

Product number	86.1252.025
----------------	-------------

Product descriptions

Name	Serological pipette, with tip, plugged, 2 ml, sterile, 25 piece(s)/bag
Product number	86.1252.025
Product description	Serological pipette, with tip, plugged, 2 ml, coding colour: green, material: PS, sterile, 25 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Diameter	6 mm
Length of product	270 mm

Product characteristics

Design	with tip, plugged
Color of print/label	green
Graduation	yes

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	25
Piece(s) / case	1000
Piece(s) / pallet	81000
Height of case	195 mm
Width of case	250 mm
Depth of case	413 mm
Volume of case	0.0201 cbm
Weight of case	4.21 kg
EAN inner box	4038917127219
EAN case	4038917021098

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT