




Product specification

**Pipette tip, 20 µl, transparent, PCR Performance Tested,
Low retention, 96 piece(s)/box**



Product information

Product number	70.3020.210
----------------	-------------

Product descriptions

Name	Pipette tip, 20 µl, transparent, PCR Performance Tested, Low retention, 96 piece(s)/box
Product number	70.3020.210
Product description	Pipette tip, working volume: 20 µl, transparent, filling level rings, PCR Performance Tested, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette and Brand and products of identical design, 96 piece(s)/box

Material & colours

Color of product	transparent
------------------	-------------

Size

Length of product	46 mm
Max. volume	20 µl

Product characteristics

Filter	no
--------	----

Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	96000
Height of inner box	148 mm
Width of inner box	111 mm
Depth of inner box	330 mm
Height of case	245 mm
Width of case	325 mm
Depth of case	348 mm
Volume of case	0.0277 cbm
Weight of case	2.82 kg
EAN inner box	4038917894197
EAN case	4038917894203

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT