



Product specification

Sarpette® M, single-channel, variable, 0.5 µl - 10 µl, grey, 1 piece(s)/case



Product information

Product number	90.3100.010
----------------	-------------

Product descriptions

Name	Sarpette® M, single-channel, variable, 0.5 µl - 10 µl, grey, 1 piece(s)/case
Product number	90.3100.010
Product description	Sarpette® M, single-channel, variable, 0.5 µl - 10 µl, grey, 1 piece(s)/case

Material & colours

Color of product	grey
------------------	------

Product characteristics

Autoclavable (max. °C)	121 °C
------------------------	--------

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Height of inner box	30 mm
---------------------	-------

Width of inner box	90 mm
Depth of inner box	315 mm
Height of case	30 mm
Width of case	90 mm
Depth of case	320 mm
Volume of case	0.0009 cbm
Weight of case	0.18 kg
EAN inner box	4038917936729
EAN case	4038917936736

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT