



Product specification

**Sarpette® M, single-channel, variable, 20 µl - 200 µl,
yellow, 1 piece(s)/case**



Product information

Product number	90.3100.200
----------------	-------------

Product descriptions

Name	Sarpette® M, single-channel, variable, 20 µl - 200 µl, yellow, 1 piece(s)/case
Product number	90.3100.200
Product description	Sarpette® M, single-channel, variable, 20 µl - 200 µl, yellow, 1 piece(s)/case

Material & colours

Color of product	grey
------------------	------

Product characteristics

Autoclavable (max. °C)	121 °C
------------------------	--------

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Height of inner box	30 mm
Width of inner box	90 mm

Depth of inner box	315 mm
EAN inner box	4038917937009
EAN case	4038917937016

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT