



## Product specification

**Sarpette® M, single-channel, variable, 20 µl - 200 µl,  
yellow, 1 piece(s)/case**



## Product information

|                |             |
|----------------|-------------|
| Product number | 90.3100.200 |
|----------------|-------------|

## Product descriptions

|                     |   |
|---------------------|---|
| Name                | Sarpette® M, single-channel, variable,<br>20 µl - 200 µl, yellow, 1 piece(s)/case |
| Product number      | 90.3100.200   |
| Product description | Sarpette® M, single-channel, variable,<br>20 µl - 200 µl, yellow, 1 piece(s)/case |

## Material & colours

|                  |      |
|------------------|------|
| Color of product | grey |
|------------------|------|

## Product characteristics

|                        |        |
|------------------------|--------|
| Autoclavable (max. °C) | 121 °C |
|------------------------|--------|

## Purity & certification

|                  |               |
|------------------|---------------|
| Product category | not regulated |
| CE certified     | no CE mark    |
| Batched          | no            |

## Packaging

|                     |       |
|---------------------|-------|
| Height of inner box | 30 mm |
|---------------------|-------|

|                    |               |
|--------------------|---------------|
| Width of inner box | 90 mm         |
| Depth of inner box | 315 mm        |
| Height of case     | 30 mm         |
| Width of case      | 90 mm         |
| Depth of case      | 320 mm        |
| Volume of case     | 0.0009 cbm    |
| Weight of case     | 0.18 kg       |
| EAN inner box      | 4038917937009 |
| EAN case           | 4038917937016 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT