



## Product specification

**Sarpette® M, single-channel, variable, 0.5 ml - 5 ml, violet,  
1 piece(s)/case**



## Product information

|                |             |
|----------------|-------------|
| Product number | 90.3100.555 |
|----------------|-------------|

## Product descriptions

|                     |   |
|---------------------|---|
| Name                | Sarpette® M, single-channel, variable, 0.5 ml - 5 ml, violet, 1 piece(s)/case |
| Product number      | 90.3100.555   |
| Product description | Sarpette® M, single-channel, variable, 0.5 ml - 5 ml, violet, 1 piece(s)/case |

## Material & colours

|                  |      |
|------------------|------|
| Color of product | grey |
|------------------|------|

## Product characteristics

|                        |        |
|------------------------|--------|
| Autoclavable (max. °C) | 121 °C |
|------------------------|--------|

## Purity & certification

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | no                      |

## Packaging

|                     |       |
|---------------------|-------|
| Height of inner box | 45 mm |
| Width of inner box  | 88 mm |

|                    |               |
|--------------------|---------------|
| Depth of inner box | 246 mm        |
| EAN inner box      | 4038917937146 |
| EAN case           | 4038917937153 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**