



Product specification

**Sarpette® M, single-channel, variable, 0.5 ml - 5 ml, violet,
1 piece(s)/case**



Product information

Product number	90.3100.555
----------------	-------------

Product descriptions

Name	Sarpette® M, single-channel, variable, 0.5 ml - 5 ml, violet, 1 piece(s)/case
Product number	90.3100.555
Product description	Sarpette® M, single-channel, variable, 0.5 ml - 5 ml, violet, 1 piece(s)/case

Material & colours

Color of product	grey
------------------	------

Product characteristics

Autoclavable (max. °C)	121 °C
------------------------	--------

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

Packaging

Height of inner box	45 mm
---------------------	-------

Width of inner box	88 mm
Depth of inner box	246 mm
Height of case	50 mm
Width of case	90 mm
Depth of case	250 mm
Volume of case	0.0011 cbm
Weight of case	0.215 kg
EAN inner box	4038917937146
EAN case	4038917937153

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT