



## Product specification

**Sarpette® M, 8-channel, variable, 20 µl - 200 µl, yellow, 1 piece(s)/case**



## Product information

Product number	90.3108.200
----------------	-------------

## Product descriptions

Name	Sarpette® M, 8-channel, variable, 20 µl - 200 µl, yellow, 1 piece(s)/case
Product number	90.3108.200
Product description	Sarpette® M, 8-channel, variable, 20 µl - 200 µl, yellow, 1 piece(s)/case

## Material & colours

Color of product	grey
------------------	------

## Product characteristics

Autoclavable (max. °C)	121 °C
------------------------	--------

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Height of inner box	28 mm
---------------------	-------

Width of inner box	150 mm
Depth of inner box	288 mm
Height of case	28 mm
Width of case	150 mm
Depth of case	288 mm
Volume of case	0.0012 cbm
Weight of case	0.455 kg
EAN inner box	4038917937337
EAN case	4038917937344

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT