



## Product specification

**Sarpette® M, 12-channel, variable, 30 µl - 300 µl, orange, 1 piece(s)/case**



## Product information

Product number	90.3112.300
----------------	-------------

## Product descriptions

Name	Sarpette® M, 12-channel, variable, 30 µl - 300 µl, orange, 1 piece(s)/case
Product number	90.3112.300
Product description	Sarpette® M, 12-channel, variable, 30 µl - 300 µl, orange, 1 piece(s)/case

## Material & colours

Color of product	grey
------------------	------

## Product characteristics

Autoclavable (max. °C)	121 °C
------------------------	--------

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Height of inner box	28 mm
---------------------	-------

Width of inner box	150 mm
Depth of inner box	288 mm
Height of case	30 mm
Width of case	150 mm
Depth of case	290 mm
Volume of case	0.0013 cbm
Weight of case	0.52 kg
EAN inner box	4038917937740
EAN case	4038917937757

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT