



## Product specification

**Urine container NFT, 3 l, with integrated needle-free transfer unit, brown, with light protection, g**



## Product information

Product number	77.575.900
----------------	------------

## Product descriptions

Name	Urine container NFT, 3 l, with integrated needle-free transfer unit, brown, with light protection, graduated
------	--

Product number	77.575.900
----------------	------------

Product description	Urine container NFT, Urine collection with hygienically closed, needle-free collection function for Urine Monovettes, max. working volume: 3 l, Ø opening: 83.5 mm, brown, with light protection, graduated, material: PE, with paper label, german/english, cap: white, with integrated needle-free transfer unit, assembled, 30 piece(s)/case
---------------------	---

## Material & colours

Color of product	brown
------------------	-------

Closure material	Polypropylene (PP)
------------------	--------------------

Color of cap	white
--------------	-------

## Size

Length excluding cap	243 mm
----------------------	--------

Length including cap	246 mm
Width of product	130 mm
Length of product	125 mm
Diameter of opening	83.5 mm
Max. volume	3 l
Sample volume	3000 ml

## Product characteristics

Design	with integrated needle-free transfer unit
Label/Print	With paper label
Color of print/label	white/black
Graduation	yes
Light protection	yes
Viewing window	yes
Cap	Membrane cap
Closure type	Screw cap
Connection	Luer (m)
Example of use	closed urine transfer into a Monovette® Urine

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes
Safety data sheet	no

## Packaging

Piece(s) / inner box	30
Piece(s) / case	30
Piece(s) / pallet	270
Height of inner box	515 mm
Width of inner box	403 mm
Depth of inner box	643 mm
Height of case	403 mm
Width of case	525 mm
Depth of case	643 mm
Volume of case	0.136 cbm
Weight of case	5.1 kg
EAN inner box	4038917947718
EAN case	4038917947725

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT