
×

Product specification

Filter tip, 100 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/box



Product information

Product descriptions

Name	Filter tip, 100 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/Box
Product number	70.3030.275
Product description	Filter tip, working volume: 100 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/box

Material & colours

Color of product	transparent
Filter material	Polyethylene (PE)

Size

Length of product	51 mm
Max. volume	100 μΙ



Filter	yes
Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Batched	yes

Packaging

480
1920
96000
148 mm
111 mm
330 mm
245 mm
325 mm
348 mm
0.0277 cbm
5.161 kg
4038917846165
4038917846172

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

