




## Product specification

**Filter tip, 200 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/SingleRefill**



## Product information

Product number	70.3031.375
----------------	-------------

## Product descriptions

Name	Filter tip, 200 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/SingleRefill
Product number	70.3031.375
Product description	Filter tip, working volume: 200 µl, transparent, filling level rings, Biosphere® plus, Low retention, suitable for SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit and Brand and products of identical design, 96 piece(s)/SingleRefill

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Filter material	Polyethylene (PE)

## Size

Length of product	51 mm
Max. volume	200 µl

## Product characteristics

Filter	yes
Graduation	yes

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Batched	yes

## Packaging

Minimum order quantity	1920
Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	115200
Height of inner box	148 mm
Width of inner box	105 mm
Depth of inner box	291 mm
Height of case	170 mm
Width of case	307 mm
Depth of case	437 mm
Volume of case	0.0228 cbm
Weight of case	5.161 kg
EAN inner box	4038917845885
EAN case	4038917845892

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT