



Product specification

SafeSeal reaction tube, 5 ml, PP



Product information

| | |
|----------------|------------|
| Product number | 72.701.500 |
|----------------|------------|

Product descriptions

| | |
|---------------------|---|
| Name | SafeSeal reaction tube, 5 ml, PP |
| Product number | 72.701.500 |
| Product description | SafeSeal reaction tube, working volume: 5 ml, material: PP, transparent, cap: natural, SafeSeal cap, cap attached, with injected graduation and writing space, 100 piece(s)/bag |

Material & colours

| | |
|------------------|--------------------|
| Color of product | transparent |
| Closure material | Polypropylene (PP) |
| Color of cap | natural |

Size

| | |
|-------------------|-------|
| Diameter | 17 mm |
| Length of product | 60 mm |

Product characteristics

| | |
|--------------------------|--------------|
| Graduation | yes |
| Cap | SafeSeal cap |
| Centrifugation max (RCF) | 25000 x g |

Purity & certification

| | |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 200 |
| Piece(s) / case | 800 |
| Piece(s) / pallet | 38400 |
| Height of inner box | 124 mm |
| Width of inner box | 184 mm |
| Depth of inner box | 294 mm |
| Height of case | 145 mm |
| Width of case | 398 mm |
| Depth of case | 598 mm |
| Volume of case | 0.0345 cbm |
| Weight of case | 3.5104 kg |
| EAN inner box | 4038917952514 |
| EAN case | 4038917952521 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT