



## Product specification

### Pressure filtration unit, PES, pore size: 0.2 µm, for sterile filtration



## Product information

|                |         |
|----------------|---------|
| Product number | 83.3944 |
|----------------|---------|

## Product descriptions

|                     |   |
|---------------------|---|
| Name                | Pressure filtration unit, PES, pore size: 0.2 µm, for sterile filtration  |
| Product number      | 83.3944   |
| Product description | Pressure filtration unit, Filtropur L 0.2 S, pre-filter: GF, membrane: PES, filtration surface: 20 cm², pore size: 0.2 µm, for sterile filtration, connection: Tube connector, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag |

## Material & colours

|                   |                        |
|-------------------|------------------------|
| Membrane material | Polyethersulfone (PES) |
|-------------------|------------------------|

## Size

|                   |          |
|-------------------|----------|
| Filtration volume | 10000 ml |
|-------------------|----------|

## Product characteristics

|                    |        |
|--------------------|--------|
| Filter             | yes    |
| Area of the filter | 20 cm² |
| Pore size          | 0.2 µm |

|                |                        |
|----------------|------------------------|
| Connection     | Tube connector         |
| Example of use | for sterile filtration |

## Purity & certification

|                  |  |
|------------------|--|
| Product category | not regulated  |
| CE certified     | no CE mark   |
| Purity standard  | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation    | Gas sterilisation                                    |
| Batched          | yes  |

## Packaging

|                      |               |
|----------------------|---------------|
| Piece(s) / inner box | 10            |
| Piece(s) / case      | 50            |
| Piece(s) / pallet    | 400           |
| Height of inner box  | 123 mm        |
| Width of inner box   | 155 mm        |
| Depth of inner box   | 352 mm        |
| Height of case       | 280 mm        |
| Width of case        | 370 mm        |
| Depth of case        | 530 mm        |
| Volume of case       | 0.0549 cbm    |
| Weight of case       | 4.24 kg       |
| EAN inner box        | 4038917926225 |
| EAN case             | 4038917926232 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT