



Product specification

Pressure filtration unit, PES, pore size: 0.2 µm, for sterile filtration



Product information

Product number	83.3944.001
----------------	-------------

Product descriptions

Name	Pressure filtration unit, PES, pore size: 0.2 µm, for sterile filtration
Product number	83.3944.001
Product description	Pressure filtration unit, Filtropur L 0.2 LS, pre-filter: GF, membrane: PES, filtration surface: 20 cm², pore size: 0.2 µm, for sterile filtration, connection: Luer lock, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Material & colours

Product material	Polycarbonate (PC)
Membrane material	Polyethersulfone (PES)

Size

Filtration volume	10000 ml
-------------------	----------

Product characteristics

Filter	yes
Area of the filter	20 cm²

Pore size	0.2 µm
Connection	Luer lock
Example of use	for sterile filtration

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Minimum order quantity	10
Piece(s) / inner box	10
Piece(s) / case	50
Piece(s) / pallet	400
Height of inner box	123 mm
Width of inner box	155 mm
Depth of inner box	352 mm
Height of case	280 mm
Width of case	370 mm
Depth of case	530 mm
Volume of case	0.0549 cbm
Weight of case	4.24 kg
EAN inner box	4038917926188
EAN case	4038917926195

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT