




Product specification

Urine transfer unit for V-Monovette® Urine, 50 piece(s)/minigrip bag



Product information

Product number	11.1240.100
----------------	-------------

Product descriptions

Name	Urine transfer unit for V-Monovette® Urine, 50 piece(s)/minigrip bag
Product number	11.1240.100
Product description	Urine transfer unit for V-Monovette® Urine, for Urine transfer, 50 piece(s)/minigrip bag

Material & colours

Color of product	yellow
------------------	--------

Size

Length of tube	10 cm
----------------	-------

Product characteristics

Design	With wings on the tube holder
Light protection	no
Connection	Vacuum tube
Type of collection	Vacuum

Purity & certification

Product category	diagnostic device, CE
CE certified	CE supplier
Batched	yes
Safety data sheet	no

Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	30000
Height of inner box	10 mm
Width of inner box	250 mm
Depth of inner box	300 mm
Height of case	170 mm
Width of case	350 mm
Depth of case	380 mm
Volume of case	0.0226 cbm
Weight of case	1.8 kg
EAN inner box	6947462401157
EAN case	6947462401423

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT